

FOCUSED ION BEAM SYSTEM: HITACHI FB-2100

TEM SAMPLE PREPARATION PROCEDURE

This procedure assumes the following:

1. Bulk sample is loaded on the SEM stage.
2. **HV** is on and the tip is conditioned

Procedure

1. Find region of interest:
 - a. Select viewing beam 40-0-30.
 - b. Use stage controls to pan to target area.
 - c. Adjust stage Z height to set focus. If the sample point is focused on the 40-0-30 beam, then all the beams will be in focus.
2. Click **Fabrication** button at the top of the software interface to open the fabrication sub-program.
3. Set **magnification** to 2.0K.
4. Click the square icon at the top-right corner of the left-hand panel. This captures an image on which to overlay the fabrication diagram.
5. Select your deposition and sputtering regions:
 - d. Use **File->Open** to open a pre-saved cutting scheme (fib-ws), or
 - e. Use the **Deposition** and **Sputtering** buttons at the top of the **Fabrication** sub-program to draw boxes that represent regions to be cut or deposited upon.

For details of the recommended deposition and sputtering regions to prepare a TEM sample, refer to the document titled "FB-2100 Micro Sampling High-Speed Processing Procedure".

When you create a "sputtering figure" or "deposition figure" you must configure a number of parameters:

- i. Interlace: set to "none" for sputtering, and "8" for deposition.
- ii. Dwell time: set to 3.0 for sputtering, and 0.5 for deposition. It is also possible to set it higher for sputtering (e.g. 75) and do a single pass (see v. below)
- iii. Width/height: as needed
- iv. Material: as needed

- v. Time/Depth/N: determines how deep the cut is made. If your “material” is available as an option, you can use the “depth” setting. Otherwise you must select a cutting time or number of passes (N).
 - vi. Scan: set scan direction and pattern.
- f. Use **Tools->Pattern List** to re-arrange the cutting order of sputtering and deposition patterns, as needed.
- g. Select the figures you want to include in the fabrication run. Click the “double arrow” icon at the top of the fabrication panel and:
- i. Drag a box around the figures, or
 - ii. Left-click the figures, one by one
6. Click the blue **fabricate** icon to begin the production run.